

SPECIFICATIONS:

- The specimens to be ground in 60, 80 or 120 grit after cutting operation
- Belt Grinder, for flat, course grinding prior to fine grinding for Metallographic Specimens
- Table top model covered with Powder Coated sheets.
- 1/2 HP Single Phase Motor
- 100mm x 915mm Endless Belt
- Easy Belt replacement and built in centering mechanism.

MODELS:

- BS1
- WBS1



- For each specimen mounting, 10-15 minutes is enough and a mirror-like specimen can be acquired.
- For permanent mounting of metallographic or rock specimen in thermostatic Plastic.
- The unit is of a rugged construction and is activated by a manual hydraulic press, capacity 6-8 ton.
- The Press is supplied complete with mould heater and sensor.
- As a standard outfit 25 or 38 mm dia mould assembly is supplied along with the machine.
- The temperature is controlled with digital controller.
- Timer facility is given within the panel

MODELS:

MP1
MP2



Single Disc Polishing Machine

SPECIFICATIONS:

- Independent Motor for Each Disc Powered by High Torque 1/2hp Motor.
- It has Water Faucet and STAINLESS STEEL Sink with Smooth Drainage System.
- It has VARIABLE SPEED SYSTEM operated manually by handle.
- It has GUN METAL DISC of Dia 200mm with LOCKING RING TO Hold the Cloth or Emery Paper.
- It is Covered with POWER COATED Sheets.

MODELS:

• SDP1 • DDP2

Note: As we are continuously improving & developing our product, the equipment supplied may be advanced in details with the description and/or illustrations shown in this catalogue.



Double Disc Polishing Machine





BELT SANDER

INVERTED METALLURGICAL

MICROSCOPE



SINGLE DISC

POLISHING MACHINE

PROFILE PROJECTOR

SPECIMEN MOUNTING PRESS

We Set **Complete Quality Metal Testing Lab**

Office: 257, Pocket A-4, Sector-4, Rohini, Delhi-110085,

Works: RZ-23/C, Vishnu Garden, Gali No.-4, Delhi-110018. Mob: 09212429102, 09810375025, Ph: 011-27044606

E-mail: vishalopticalworks@gmail.com, Web: www.vishalopticalworks.com